

Key improvements of the Binder Lift 2.0



72% smaller storage footprint*



42% lighter*



Softer and more pliable material

Model Comparison

Nylon

(Discontinued)



Vinyl

(Discontinued)



2.0



Weight			
Standard Size	4lbs	5.5lbs	3.2lbs
Bariatric Size	4.5lbs	7lbs	3.6lbs
Storage Dimensions			
Standard Size	16"x13"x4"	24"x16"x5"	12.5"x9.5"x4.5"
Bariatric Size	16"x13"x5"	24"x17"x5"	13.5"x9.5"x6"
Material			
Torso Material	Denier 1,000 Nylon	18oz Vinyl	14oz Vinyl
Handle Material	2" polypropylene	1.5" TPU Coated Webbing	1.5" PVC Coated Webbing
Leg Strap Material	2" TPU Coated Webbing	2" TPU Coated Webbing	2" PVC Coated Webbing

*compared to Vinyl model